

# PBAIR®

## IAPV ACCESSORIES

### (INTERMITTENT ABDOMINAL PRESSURE VENTILATION)

The PBAir® is the new and innovative interface that allows the Intermittent Abdominal Pressure Ventilation.


The use of PBAir® with a dedicated ventilator software (BB respirator) allows the lung ventilation (increasing gas-exchange) in patients with muscular weakness and diaphragmatic paralysis, conditions of all Neuromuscular patients:



The IAPV ventilation promotes an extra-thoracic Non Invasive Ventilation, without the aid of masks, avoiding the tracheotomy, reducing the facial decubitus risk due to the mask, improving the patient's quality of life.

Below, sizes available for PBAir®.

NOTE ASIDE: *It is possible to provide, on request, a size that is not present in the table below.*



| Size        |        |                     |
|-------------|--------|---------------------|
| Code        | Length | Waist circumference |
| PBAir0090S  | 90 cm  | 70 cm               |
| PBAir00100S | 100 cm | 80 cm               |
| PBAir00115M | 115 cm | 95 cm               |
| PBAir00120M | 120 cm | 100 cm              |
| PBAir00128M | 128 cm | 108 cm              |
| PBAir00134M | 134 cm | 114 cm              |
| PBAir00144M | 144 cm | 124 cm              |
| PBAir00150L | 150 cm | 130 cm              |